

ACE CALIBRATOR

REF 30550

6 x1 mL

Pour d'autres langues
Für andere Sprachen
Para otras lenguas
Per le altre lingue
Dla innych języków

Para outras línguas
Για τις άλλεςλώσσες
För andra språk
For andre språk
For andre sprog



www.trinitybiotech.com

ASSIGNED VALUES

The assigned values are lot specific and can be found on the enclosed component list.

INTENDED USE

ACE Calibrator is intended to establish the point of reference in the kinetic determination of angiotensin converting enzyme (ACE) in human serum or plasma.

REAGENT

ACE Calibrator, Catalog No. 30550
Lyophilized preparation containing ACE (porcine) in a buffered human serum base with stabilizers and preservatives.

PRECAUTIONS

ACE Calibrator is for "in vitro Diagnostic Use". Normal precautions exercised in handling laboratory reagents should be followed. Dispose of waste observing all local, state and federal laws.

ACE Calibrator is a POTENTIALLY BIOHAZARDOUS MATERIAL. Source materials from which this product was derived were found negative for HBsAg and for antibodies against HCV, HIV-1 and HIV-2 by approved test methods. However, since no test method can offer complete assurance that infectious agents are absent, this product should be handled observing the same safety precautions employed when handling any potentially infectious material.

Angiotensin converting enzyme is classified as a reproductive hazard in California Proposition 65. Refer to Material Safety Data Sheet for any updated risk, hazard or safety information.

MATERIALS REQUIRED BUT NOT PROVIDED

- Deionized water
- Pipet calibrated to deliver 1.0 mL

PREPARATION

1. Add 1.0 mL water to vial of ACE Calibrator.
2. Swirl gently and allow to stand for 5 minutes.
3. Invert gently and mix well until dissolution is complete. Swirl gently to mix before each use.

STORAGE AND STABILITY

Store ACE Calibrator in the refrigerator (2–8°C). Vial label bears expiration date.

After reconstitution, ACE Calibrator solution is stable for at least 7 days in refrigerator (2–8°C).

DETERIORATION

Do not use ACE Calibrator if there is evidence of microbial growth.










GENERAL INSTRUCTIONS FOR USE

The ACE Calibrator must be included with test samples in the determination of angiotensin converting enzyme by Trinity Biotech Procedure No. 305-UV. Refer to the Trinity Biotech procedure for its use in calculating the ACE activity of tests.

Trinity Biotech, Inc., warrants that its products conform to the information contained in this and other Trinity Biotech publications. Purchaser must determine the suitability of the product for its particular use.

ORDERING INFORMATION

KIT	Item	ACE Calibrator
Catalog No.	Item	Quantity
30550	ACE Calibrator	6x1 mL
	ACE Controls	
A6040	Normal	6x1 mL
A7040	Elevated	6x1 mL

 Manufactured	CONTROL N Normal Control	CONT Contents
EC REP Authorized Representative	CONTROL I Intermediate Control	SUBS Substrate
 Consult accompanying documents	CONTROL E Elevated Control	REAG Reagent
REF Product Number	CAL Calibrator	DIL Diluent
LOT Lot	C.F. Coefficient Factor	ACT Activator
 Use by	RNG Range	IVD For In Vitro Diagnostic use
 Caution, consult accompanying documents	STD Standard	 Biohazard
 Store at 2–8°C	 Store at 2–30°C	 or Hazard
 Trinity Biotech USA Jamestown, NY 14701 Tel. 1 800-325-3424 Fax: 716-488-1990	EC REP CE	Trinity Biotech plc Bray Co. Wicklow, Ireland Tel. 353 1 2769800 Fax 353 1 2769888 www.trinitybiotech.com

30550-29 Rev E
10/2011